

PurSolv Ultra Plus


High-Purity Solvents for PFAS Testing

Why PFAS Testing Matters

Per- and polyfluoroalkyl substances (PFAS) are synthetic chemicals known for their persistence, toxicity, and environmental impact. Accurate trace-level detection of PFAS demands solvents of exceptional purity - particularly methanol, acetonitrile, and water - used in LC-MS/MS analysis.

Our PFAS-tested solvents are engineered to eliminate background interference and prevent false positives, ensuring reliable and compliant PFAS testing.

PFAS Structure



Product Features:

- PFAS-Free certification as per EPA Methods 533 and 537.1
- Lot-to-Lot consistency
- Interference-free baseline for LC-MS/MS chromatograms
- Low background in both ESI/APCI (+) and (-) ionization modes
- Controlled metal ion levels (<5 ppb)
- Microfiltration & packaging in borosilicate glass bottles

Applications:

Drinking water
testing (EPA 533
& 537.1)

Wastewater,
groundwater &
surface water
analysis

Landfill leachate
and food chain
monitoring

Method validation
and instrument
calibration

Product Range

Product No.	Description	Cas No.	Packing
PSR49378	Acetonitrile PurSolv Ultra Plus tested for EPA533 and EPA537.1 PFAS	75-05-8	1L
PSR49354	Methanol PurSolv Ultra Plus tested for EPA533 and EPA537.1 PFAS	67-56-1	1L
PSR49283	Water PurSolv Ultra Plus tested for EPA533 and EPA537.1 PFAS	7732-18-5	1L

Quality You Can Trust

Our PFAS-tested solvent range is designed to deliver consistent purity, analytical accuracy, and environmental responsibility, enabling laboratories to meet regulatory expectations with total confidence.

PureSynth Research Chemicals GmbH

64683 Marie-Curie-Str. 3, Einhausen, Germany

☎ 1800-8908-260 ✉ sales@pure-synth.com 🌐 www.pure-synth.com